

# REMEDIC

## Quick Fixing Bandage

**re-usable, self-adhesive protective bandage**



- REMEDIC Quick Fixing bandages are made of 100% polyester.
- A special weave allows the lighter underside stick to the darker top side.
- The adhesive strength increases with the overlapping area of the material
- The strong adhesion enables you to fix e.g. a hot water bottle or an ice pack to your body without losing freedom of movement.
- REMEDIC Quick Fixing bandages can serve as a compression bandage when wrapped tightly.
- The soft material is ideal for wrapping sensitive parts of the body that need to be protected gently.
- The material keeps you warm and protects against drafts, but is breathable and pleasant on the skin
- The bandage can be easily opened without any force or pressure and the layers be separated without making any noise.
- REMEDIC Quick Fixing bandages are versatile and also suitable e.g. also for wrapping sensitive objects.
- The bandages are reusable and can be used many times – they are machine washable at temperatures up to 60 ° Celsius.



## REMEDIC Quick Fixing Bandage


the self-adhesive protective bandages are available in 4 shapes



Neck / head- 100 x 15 cm




Body - 132 x 25 cm



Triangle

Side lengths 96 / 96 / 136cm



Leg- 77 x 55 cm



For additional help/info  
email:shop@remedic.co

MADE IN CHINA